



HIGH PERFORMANCE DIGITAL IMAGING
made easy

RETIGA EXi *FAST1394*


Very High Sensitivity IEEE 1394 FireWire™ Digital CCD Camera – Monochrome or Color

The QImaging Retiga EXi CCD digital camera features enhanced visible and IR quantum efficiency resulting in very high sensitivity that is ideal for demanding low light and fluorescence imaging applications. A progressive scan interline CCD sensor gives a resolution of 1.4 million pixels in a 12-bit digital output. High-speed low noise electronics provide linear digital data for rapid image capture. The IEEE 1394 FireWire™ digital interface allows ease of use and installation with a single wire requiring no framegrabber or external power supply. The Retiga EXi includes QCapture Software for for Microsoft Windows® and Mac® OS based software systems for real time image preview and capture. A **Software Development Kit (SDK)** is available for interfacing with custom software.



Retiga EXi (non-cooled)

Retiga EXi cooled



Note: Lenses are shown for illustration only and are not included.

CAMERA MODELS

Includes: IEEE 1394 FireWire™ cable, IEEE 1394 PCI card, QCapture software and access to SDK

- **Monochrome Retiga EXi Cooled** Model: RET-EXi-F-M-12-C
- **Monochrome Retiga EXi Non-cooled** Model: RET-EXi-F-M-12 CCD Digital Camera, 12-bit
- **Color Retiga EXi Cooled** Model: RET-EXi-F-CLR-12-C
- **Color Retiga EXi Non-cooled** Model: RET-EXi-F-CLR-12 CCD Digital Camera, 12-bit

CAMERA OPTIONS

- Removable **IR cutoff filter**
- **RGB Color Filter** for monochrome cameras (F-mount interface required) Refer to spec sheet for more details
- **Extended Warranty**



FEATURES

- High Quantum Efficiency
- High Resolution 1.4 Million pixel sensor
- High Speed Readout
- Low Noise Electronics
- Optional/Removable IR cutoff filter
- Flexible Exposure Control from 10µs to 17.9min
- External Sync and Trigger

BENEFITS

- Very high sensitivity for demanding low-light & fluorescent imaging
- Highly detailed, sharp images
- Previewing & focusing in real time
- 110fps in 8x8 binning & ROI
- 10fps full resolution @ 12-bits
- Ideal for automated imaging applications
- Quantitation & imaging of low light levels
- Highly focused visible range images with IR filter in place, and removable for IR applications
- Optimal Integration over a wide range of light levels
- Tight synchronization with flashlamps, automated filters, shutters & microscope stages
- Minimizes thermal noise during low light, long exposure imaging
- Increased sensitivity for quantitation & imaging of very low light levels
- Increased frame rate
- High performance imaging outside the visible range
- Simple connectivity
- Ease of use & installation
- Portability with laptop computer
- Simultaneous use of multiple cameras through a single port
- Single cable operation, no external power supply or control unit
- Choose from a large selection of life science & industrial software for microscopy, machine vision, and video streaming functions

Peltier Cooling

Binning

Extended IR Sensitivity

IEEE 1394 FireWire™ QImaging Fast 1394 Technology


Extensive application software support

RETIGA EXi FAST1394 SPECIFICATIONS

APPLICATIONS

- Brightfield, Phase Contrast and Dark-field Microscopy
- Fluorescence Microscopy
- Live Cell Imaging
- Pathology, Histology, Cytology
- Green Fluorescent Protein (GFP) Application
- FISH
- Ca⁺⁺ Ratio Analysis
- Motility and Motion Analysis
- DNA Analysis
- Metallurgical Microscopy
- Semiconductor Inspection
- Manufacturing Quality Control
- Failure Analysis
- Forensic Analysis

SPECTRAL RESPONSE



CCD SENSOR

Light Sensitive Pixels	1.4 million; 1392 x 1040
Binning Modes	2x2, 4x4, 8x8
ROI (Region Of Interest)	From 1x1 pixels up to full resolution, continuously variable in single pixel increments
Exposure/Integration Control	10µs to 17.9min in 1µs increments
Sensor Type	Sony ICX285 Progressive Scan Interline CCD, Monochrome or Colour
Pixel Size	6.45µm x 6.45µm
Linear Full Well	18,000e ⁻ ; 22,000e ⁻ in 2x2 binning
Read noise	8e ⁻
Dark Current	0.15e ⁻ /pix/s cooled
Cooling Available	Yes (optional)
Cooling Type	Peltier thermoelectric cooling to 25 degrees Celsius below ambient
Digital Output	12-bit
Readout Frequency	20, 10, 5, 2.5MHz
Frame Rate	10fps full resolution @ 12-bits, 205*fps maximum with binning and ROI functions

CAMERA

Computer Platform/Operating System	Microsoft Windows® & Mac® OS**
Digital Interface	IEEE 1394 FireWire™
Sustained Image Data Rate	40MB/s
External Trigger	TTL Input (optically coupled)
Trigger Types	Internal, Software, External
External Sync	TTL Output (optically coupled)
Gain Control	0.7 to 30 times
Offset Control	-2048 to 2047
Optical Interface	2/3", C-Mount optical format
Threadmount	1/4" – 20 Mount
Power Requirements	7 watts non-cooled; 13 watts cooled; 8-24V
Weight	640g non-cooled; cooled 920g
Warranty	2 years
Operating environment	0 to 50 degrees Celsius (32 to 122F); Storage Temperature: -10 to 60 degrees Celsius
Humidity	Less than 80% non-condensing at 35 degrees Celsius (95F)

*Special Order Only. Standard model achieves 165fps.

**Refer to QImaging website for detailed listing of supported operating systems.

Note: Specifications are nominal and subject to change.

04-0002C-B



GT Vision Ltd
Cherry Gdns Ind Est, Helions Bumpstead Rd, Haverhil, Suffolk, CB9 7AA, UK
Tel +44(0)1440 714737
Fax +44(0)1440 709421
sales@gt-vision.com

WWW.GT-VISION.COM