

BA200 Motic Biological Microscopes from GX Optical



BA200 Binocular

PA23B49201

- Binocular head Siedentopf type 30° inclined, 360° rotating
- Widefield high eyepoint eyepiece WF10X / 18mm
- Reversed quadruple nosepiece
- CCIS E Plan Achromat objectives EPL4X, 10X, 40X S, 100X S-oil
- Coaxial coarse and fine focusing adjustment
- Built in low position coaxial mechanical stage
- Focusable 1.25 N.A. Abbe condenser
- Iris diaphragm with filter holder
- Built in halogen illumination 6V/20W with intensity control
- Universal power supply 100-240V, CE



BA200 Trinocular

PA02B49301

- Trinocular head Siedentopf type 30° inclined, 360° rotating
- Widefield high eyepoint eyepiece WF10x / 18mm
- Reversed quadruple nosepiece
- CCIS E Plan Achromat objectives EPL4X, 10X, 40X S, 100X S-oil
- Coaxial coarse and fine focusing adjustment
- Built in low position coaxial mechanical stage
- Focusable 1.25 N.A. Abbe condenser
- Iris diaphragm with filter holder
- Built in halogen illumination 6V/20W with intensity control
- Universal Power supply 100-240V, CE

Basic Accessories

EYEPIECE LENSES (per piece)

SG01S0105	WF10X / 18mm
SG02S0205	WF15X / 12mm
SG02S0405	WF20X / 11mm

PHOTO AND VIDEO ADAPTERS

SP100350	CCD adapter 1X
SP100384	CCD adapter 0.65X
SP100294	Photo adapter (requires one of the photo eyepiece below)
SG02S1001	Photo eyepiece 2.5X
SG02S1101	Photo eyepiece 4X

OBJECTIVE LENSES

SG01S012D4	CCIS E Plan Achromat 4x / 0.1
SG01S022D4	CCIS E Plan Achromat 10x / 0.25
SG01S042D4	CCIS E Plan Achromat 40x / 0.65 / S
SG01S062D4	CCIS E Plan Achromat 100x / 1.25 / S-Oil
SG01S012B5	CCIS EF Plan Achromat 4x / 0.1
SG01S022B5	CCIS EF Plan Achromat 10x / 0.25
SG01S042B5	CCIS EF Plan Achromat 40x / 0.65 / S
SG01S052B5	CCIS EF Plan Achromat 60X / 0.80 / S
SG01S062B5	CCIS EF Plan Achromat 100x / 1.25 / S-Oil

PHASE CONTRAST

SW01238A	Simple phase contrast attachment 10x (centerable)
SG01S022C1	Phase contrast objective E Plan Ph10x
SW01238B	Simple phase contrast attachment 40x (centerable)
SW0123AQ	Phase contrast objective E Plan Ph40x
SP990001	Centering Telescope

DARKFIELD

SW012391B	Simple darkfield attachment (10x - 40x)
-----------	---

SIMPLE POLARIZING

SW0199UK	Simple add on analyzer
SW0199UH	Simple add on polarizer

OPTIONAL ACCESSORIES

SG060728	Yellow filter (45mm diameter)
SG060726	Green filter (45mm diameter)
SG060727	Blue filter (45mm diameter)
SG060768	Neutral density filter N.D.2 (45mm diameter)

REPLACEMENTS

SE060096	Halogen Lamp 6V/20W
MF02020	Immersion Oil 5ml

GX Optical

A divn of GT Vision Ltd
Hazel Stub Depot
Camps Road
Haverhill
Suffolk
CB9 9AF
UK
Tel: 01440 714737
Fax: 01440 709421
E-mail: sales@gxoptical.com
www.gxoptical.com