

## K Series Motic Quad Magnification Stereo Microscopes from GX Optical



### **Model K-400P**

[PK40.011.001](#)

Binocular head  
Widefield eyepieces WF10X / 23mm  
Standard objective 1.0X  
Magnification changer 6X, 12X, 25X and 50X  
Streamline plane base stand



### **Model K-400L**

[FK40.000.201](#)

Binocular head  
Widefield eyepieces WF10X / 23mm  
Standard Objective 1.0X  
Magnification changer 6X, 12X, 25X and 50X  
Streamline plane base stand with incident and transmitted halogen lights 12V / 10W with intensity control  
Main supply 220-240V

### **Model K-401P (with plan objective 1.0X)**

[PK40.611.001](#)

Binocular head  
Widefield eyepieces WF10X / 23mm  
Plan Objective 1.0X  
Magnification changer 6X, 12X, 25X and 50X  
Streamline plane base stand

### **Model K-401L (with plan objective 1.0)**

[PK40.600.101](#)

Binocular head  
Widefield eyepieces WF10X / 23mm  
Plan Objective 1.0X  
Magnification changer 6X, 12X, 25X and 50X  
Streamline plane base stand with incident and transmitted halogen light 12V / 10W with intensity control  
Main supply 220V-240V (CE)



#### **Model K-500P**

[FK50.111.001](#)

Binocular head  
Widefield eyepieces WF10X / 23mm  
Standard objective 1.0X  
Magnification changer 6.4X, 10X, 16X, 25X and 40X  
With streamline plane base stand



#### **Model K-500L**

[FK50.100.201](#)

Binocular head  
Widefield eyepieces WF10X / 23mm  
Standard Objective 1.0X  
Magnification changer 6.4X, 10X, 16X, 25X and 40X  
Streamline plane base stand with incident and transmitted halogen light 12V / 10W with intensity control  
Main supply 220V-240V (CE)



#### **Model K-700P**

[FK70.111.001](#)

Zoom binocular head  
Standard objective 1.0x  
Zoom ratio 1:5.2  
With streamline plane base stand



#### **Model K-700L**

[FK70.100.201](#)

Zoom binocular head  
Standard objective 1.0x  
Zoom ratio 1:5.2  
Incident halogen light 12V/15W and transmitted halogen light 12V/10W, with intensity control  
Main supply 220V-240V

**Heads**

11.410.101	Stereo microscope head K-400 (without incident light)
11.410.102	Stereo microscope head K-400 (with light)
11.410.103	Stereo microscope head K-400 (with bonder arm)
PK40.621.001	Stereo microscope head K-401 (without light)
PK40.622.001	Stereo microscope head K-401 (with light)
11.510.101	Stereo microscope head K-500 (without light)
11.510.102	Stereo microscope head K-500 (with light)
11.510.103	Stereo microscope head K-500 (with bonder arm)
11.710.101	Stereo microscope head K-700 (without incident light)
11.710.102	Stereo microscope head K-700 (with light)
11.710.103	Stereo microscope head K-700 (with bonder arm)

**Stands**

11.510.131	Streamline base stand (2112)
11.510.132	Streamline base stand with transmitted light 12V/10W halogen, with intensity control Main supply 220-240V (2111)
11.510.133	Universal stand (2105)
11.510.134	Universal stand, rectangular (2105S)
11.510.135	Stand 2107 articulating arm boom stand (base stand version)
11.510.136	Stand 2108 ball bearing boom stand (base stand version)
11.510.137	Stand 2109 articulating arm boom stand (table clamp version)
11.510.138	Stand 2110 ball bearing boom stand (table clamp version)

**Eyepieces (pair)**

11.500.141	Eyepiece WF 6.25x / 23 mm (1301)
11.500.145	Eyepiece WF 10x / 23 mm (2302)
11.500.143	Eyepiece WF 15x / 17.6 mm (1303)
11.500.144	Eyepiece WF 20x / 13.4 mm (1304)
11.500.146	Eyepiece WF 32x / 12 mm (1308)

**Objectives**

SG01.T10202A	Objective, 0.3x (1201A)
SG01.T01202A	Objective 0.5x (1202)
SG01.T13201A	Objective 0.625x (1205)
SG01.T11202A	Objective 1.5x (1203)
SG01.T04202A	Objective 2.0x (1204)

**Objectives (For K-700)**

SG01.T10202A	Objective 0.3x (1201A)
SG01.T01202A	Objective 0.5x (1202A)
SG01.T13201A	Objective 0.625x (1205A)
SG01.T11202A	Objective 1.5x (1203A)
SG01.T04202A	Objective 2.0x (1204A)

**Photo Equipments**

11.500.191	Photo-equipment (1501.3) for K-500 and K-700 with T2-thread, photo eyepiece
11.400.191	Photo-equipment (1501.4) for K-400 with T2-thread, photo eyepiece
11.500.192	Photoadapter with focus-adjustment (1502)

**CCD Equipments**

11.500.221	CCD-equipment (1509.3) for K-500 and K-700 with C-mount
11.400.221	CCD-equipment (1509.4) for K-400 with C-mount

**CCD & Photo Equipments**

11.500.231	CCD-equipment included, C-mount (1510.3) Combined with: Photo equipment included, T2-thread, Photo eyepiece for K-500 and K-700
11.400.231	CCD-equipment included, C-mount (1510.4) Combined with: Photo equipment included, T2-thread, Photo eyepiece for K-400

**Spare Parts**

11.500.351	Specimen-clamp
11.500.352	Insert plate, black/white 82mm diameter
11.500.353	Insert plate, frosted glass 82mm diameter
53.100.254	Dust cover (size B)
40.000.251	Dust cover 6.9 x 28 x 28cm

**Spare Bulbs**

11.510.354	Halogen bulb 12V/10 W (incident)
11.510.355	Halogen bulb 12V/10 W (transmitted)
SP12.0014	Fluorescent ring light tube 6400K

