



## SZ51

Excellent performance  
for routine stereo microscopy

---



- 5:1 zoom ratio with a zoom magnification range of 0.8x to 4.0x
- Targeted to OEM and education
- Excellent price-performance ratio

### Optics (SZ51 / SZ61)

- ComfortView eyepieces
- The best zoom ratios in their class
- Sharp, bright and high contrast optics
- Excellent image flatness with deep depth of focus
- Choice of auxiliary objectives with working distance up to 400 mm
- A choice of 45° and 60° tube inclination models for comfortable viewing

## SZ61

Top zoom performance  
with image recording

---



- Imaging and documentation capability with a trinocular head for camera output (SZ61-TR)
- Extended 6.7:1 zoom ratio with a zoom magnification range of 0.67x to 4.5x

### Other features (SZ51 / SZ61)

- Unique LED integrated reflected/transmitted illumination stand
- Wide variety of lighting systems, stands and accessories
- Zoom stopper
- ESD safety design due to special materials and surface coatings
- Choice of auxiliary objectives with working distance up to 400 mm

## SZX7

Flexible modular system  
with highest resolution

---



- ComfortView eyepieces
- Best resolving power in this class. Superior image quality is ensured up to 336x
- Outstanding 7:1 zoom ratio with a zoom magnification range of 0.8x to 5.6x
- Tilting head gives greater comfort (SZX-TBI)
- Observation methods including GFP fluorescence and oblique illumination
- Motorised focus with PC connection and control (SZX-FOA)
- Front zoom handle with zoom click feature