


GX Microscopes

product datasheet

Compound Microscope Kits and Condensers



Fluorescence kits for L2000 and L1500 series microscopes



Incident illuminator kit for L2000 series microscopes




Plan achromatic phase kit



Slider phase kit



Simple phase kit




Immersion and dry darkfield condensers

GX Microscopes supply a very wide range of attachments for our compound microscopes. The most popular kits include a fluorescence kit which includes a transformer, MV lamphouse and two high quality fluorescence objectives. Our next most popular kit is our plan achromatic phase kit which gives outstanding results with the L2000 series in particular.

In addition to the above, adding an incident illuminator, a set of metallurgical objectives and lamp control box (incident illuminator kit) to a L2000 microscope converts it into a fully functional L2020 series materials microscope.

For details of these kits and our complete range of accessories for your GX Microscope contact us direct or your local supplier.

Supplied to you by:



www.gxmicroscopes.com

Due to our policy of continuous development specifications may change without notice
June 05