



GX Microscopes

product datasheet

Eyepieces and Reticules



Supplied to you by:



Huygens Eyepieces			
Multiple	Focal distance(mm)	Exit pupil distance (mm)	Field Range (mm)
5X	50	15.2	φ18
6X	41.7	6.6	φ16
10X	25	6	φ14
15X	16.6	6.11	φ11
Eyepieces with built-in 10mm scale reticules 23.2mm			
Multiple	Focal distance(mm)	Eyepoint Height (mm)	Field range (mm)
10X	25	16.1	φ17
Widefield Eyepieces 23.2mm			
Multiple	Focal distance(mm)	Eyepoint Height (mm)	Field Range (mm)
10X	25	19.8	φ20
10X	25	19.8	φ18
10X	25	19.2	φ18
10X	25	16.1	φ18
10X(point)	25	16.1	φ18
10X	25	15.8	φ18
10X	25	16.1	φ16
10X	25	18.5	φ18
12.5X	20	16.8	φ14
15X	16.6	15	φ11
16X	15.6	14.7	φ11
20X	12.5	15.3	φ11
20X	12.5	8.8	φ11
Eyepieces with built-in 10mm scale reticules 30mm			
Multiple	Focal Distance (mm)	Eyepoint Height (mm)	Field Range (mm)
10X	25	17.9	φ17
Wide Field Eyepieces 30mm			
Multiple	Focal Distance (mm)	Eyepoint Height (mm)	Field Range (mm)
10X	25	14.9	φ20
10X	25	17.9	φ20
15X	16.6	9.3	φ15
20X	12.5	11.1	φ12
25X	10	8.3	φ9

Also available 5mm scale and crosshair reticules.

www.gxmicroscopes.com

Due to our policy of continuous development specifications may change without notice

June 05