

GX Microscopes

product datasheet

L1180 Biological Video Microscope



Specifications

		Model:	L1180
Eyepieces	Widefield	WF10X (18mm FOV)	STD
		WF16X (11mm FOV)	OPT
Objectives	Achromatic	4X/0.10	STD
		10X/0.25	STD
		40X/0.65 (sprungloaded)	STD
		100X/1.25 Oil (sprungloaded)	STD
Nosepiece	Locating on ball bearing	Quadruple	STD

The L1180 video microscope is the easiest way to combine routine microscopy with the recording and displaying of observed images.

This inexpensive microscope includes a full range of achromatic objectives and a built-in video camera ready for connection to any video monitor, VCR or image capture card in a PC.

This microscope is ideal for use in teaching and research laboratories as it provides all the features needed for routine observation combined with a very simple way to obtain video images.

Supplied to you by:

www.gxmicroscopes.com

Due to our policy of continuous development specifications may change without notice

May 04