

# GX Microscopes

## product datasheet

### L301 Biological Microscope



#### Specifications

		Model:	L301
Eyepieces	Widefield	WF10X (14mm FOV)	STD
		WF15X (11mm FOV)	STD
Objectives	Achromatic (35mm)	4X/0.10	STD
		10X/0.25	STD
		40X/0.65 (sprungloaded)	STD
		100X/1.25 Oil (sprungloaded)	STD
Total Magnification		40X – 1500X	STD
Condenser	Abbe condenser	NA=1.25 with aperture diaphragm and filter seat	STD
Illumination	Transmitted light	Concave reflecting mirror (requires light source)	STD

The L301 biological upright microscope provides great quality microscopy at an astonishingly low price.

Incorporating coarse and fine focus controls and a dual axis mechanical stage normally found on research microscopes, the user can be assured of comfortable, effortless microscopy over prolonged periods.

A complete range of achromatic objectives and an Abbe condenser complete this high specification to produce great quality images.

This robust microscope is ideal for amateur microscopes, schools and colleges and can be supplied with an optional wooden storage case.

#### Supplied to you by:

[www.gxmicroscopes.com](http://www.gxmicroscopes.com)

Due to our policy of continuous development specifications may change without notice

May 04