

# GX Microscopes

## product datasheet

### Slide Boxes

Quality, Durable Plastic  
Slide Boxes

#### Features and Specifications

GX Microscopes range of high quality, incredibly good value slide boxes. Constructed from durable, shock resistant plastic, these containers have clip-lock mechanisms to ensure a long life for your precious specimen slides.


All except the 6 slide box include a slide listing sheet.

These containers are suitable for sending slides by post but it is a good idea to reduce rattling of the slides by applying tough tape over the tops of the slides before shipment. Also check that the foam is still in a good condition before sending by post.

All slide boxes are supplied in white. Other colours are available for large volume customers.



Supplied to you by:



[www.gxmicroscopes.com](http://www.gxmicroscopes.com)

Due to our policy of continuous development specifications may change without notice

Jan 06