

# GX Microscopes

product datasheet

## Quality

Microscope Slides · Coverslips · Prepared Slides

### GXM Glass Slides Specifications

#### Microscope Slides



7101 & 7102



7103



7104



7105



7105-S







Box of 72

| Item  |  | Specification and Packing |                    |
|---|--|---------------------------|--------------------|
| <b>7101</b><br>Plain Slide<br>Ground Edges          |  | Size mm                   | 25.4 X 76.2        |
|   |  | Thickness(mm)             | 1.0-1.2<br>0.8-1.0 |
|   |  | Packing (PCS/Box)         | 50<br>72           |
|   |  | Boxes per Carton          | 50                 |
| <b>7102</b><br>Plain Slide<br>Unground Edges        |  | Size mm                   | 25.4 X 76.2        |
|   |  | Thickness(mm)             | 1.0-1.2<br>0.8-1.0 |
|   |  | Packing (PCS/Box)         | 50<br>72           |
|   |  | Boxes per Carton          | 50                 |
| <b>7103</b><br>Single Concave Slide<br>Ground Edges |  | Size mm                   | 25.4 X 76.2        |
|   |  | Thickness(mm)             | 1.1-1.3            |
|   |  | Packing (PCS/Box)         | 50<br>72           |
|   |  | Boxes per Carton          | 50                 |

# GXM Glass Slides Specifications


## Microscope Slides



| Item  |   | Specification and Packing |                    |
|---|---|---------------------------|--------------------|
| <b>7104</b><br>Double Concave Slide<br>Ground Edges                             |    | <b>Size mm</b>            | 25.4 X 76.2        |
|   |   | <b>Thickness(mm)</b>      | 1.1-1.3            |
|   |   | <b>Packing (PCS/Box)</b>  | 50<br>72           |
|   |   | <b>Boxes per Carton</b>   | 50                 |
| <b>7105</b><br>Markable Slide<br>Frosted<br>One End One Side<br>Ground Edges    |   | <b>Size mm</b>            | 25.4 X 76.2        |
|   |   | <b>Thickness(mm)</b>      | 1.0-1.2<br>0.8-1.0 |
|   |   | <b>Packing (PCS/Box)</b>  | 50<br>72           |
|   |   | <b>Boxes per Carton</b>   | 50                 |
| <b>7105-1</b><br>Markable Slide<br>Frosted<br>One End On Side<br>Unground Edges |  | <b>Size mm</b>            | 25.4 X 76.2        |
|   |   | <b>Thickness(mm)</b>      | 1.0-1.2<br>0.8-1.0 |
|   |   | <b>Packing (PCS/Box)</b>  | 50<br>72           |
|   |   | <b>Boxes per Carton</b>   | 50                 |

# GXM Glass Slides Specifications


## Microscope Slides



| Item   |  | Specification and Packing |                    |
|--|--|---------------------------|--------------------|
| <b>7106</b><br>Markable Slide<br><br>Frosted Both Ends<br><br>One Side<br><br>Ground Edges |  | <b>Size mm</b>            | 25.4 X 76.2        |
|  |  | <b>Thickness(mm)</b>      | 1.0-1.2<br>0.8-1.0 |
|  |  | <b>Packing (PCS/Box)</b>  | 50<br>72           |
|  |  | <b>Boxes per Carton</b>   | 50                 |
| <b>7107</b><br>Markable Slide<br><br>Frosted One End Both Sides<br><br>Ground Edges        |  | <b>Size mm</b>            | 25.4 X 76.2        |
|  |  | <b>Thickness(mm)</b>      | 1.0-1.2<br>0.8-1.0 |
|  |  | <b>Packing (PCS/Box)</b>  | 50<br>72           |
|  |  | <b>Boxes per Carton</b>   | 50                 |
| <b>7107-1</b><br>Markable Slide<br><br>Frosted One End Both Sides<br><br>Unground Edges    |  | <b>Size mm</b>            | 25.4 X 76.2        |
|  |  | <b>Thickness(mm)</b>      | 1.0-1.2<br>0.8-1.0 |
|  |  | <b>Packing (PCS/Box)</b>  | 50<br>72           |
|  |  | <b>Boxes per Carton</b>   | 50                 |

# GXM Glass Slides Specifications

## Microscope Slides



| Item  |  | Specification and Packing |                    |
|---|--|---------------------------|--------------------|
| <b>7108</b><br>High-Quality Slide<br><br>Frosted<br>Both Ends<br>Both Sides<br><br>Ground Edges       |  | <b>Size mm</b>            | 25.4 X 76.2        |
|   |  | <b>Thickness(mm)</b>      | 1.0-1.2<br>0.8-1.0 |
|   |  | <b>Packing (PCS/Box)</b>  | 50<br>72           |
|   |  | <b>Boxes per Carton</b>   | 50                 |
| <b>7109</b><br>High Quality Slide<br><br>Color Frosted<br><br>One End<br>One Side<br><br>Ground Edges |  | <b>Size mm</b>            | 25.4 X 76.2        |
|   |  | <b>Thickness(mm)</b>      | 1.0-1.2<br>0.8-1.0 |
|   |  | <b>Packing (PCS/Box)</b>  | 50<br>72           |
|   |  | <b>Boxes per Carton</b>   | 50                 |
| <b>7110</b><br>Markable Slide<br><br>Frosted<br>One side<br><br>Ground Edges                          |  | <b>Size mm</b>            | 25.4 X 76.2        |
|   |  | <b>Thickness(mm)</b>      | 1.0-1.2<br>0.8-1.0 |
|   |  | <b>Packing (PCS/Box)</b>  | 50<br>72           |
|   |  | <b>Boxes per Carton</b>   | 50                 |
| <b>7111</b><br>Markable Slide<br><br>Frosted<br>Both Sides<br><br>Ground Edges                        |  | <b>Size mm</b>            | 25.4 X 76.2        |
|   |  | <b>Thickness(mm)</b>      | 1.0-1.2            |
|   |  | <b>Packing (PCS/Box)</b>  | 50<br>72           |
|   |  | <b>Boxes per Carton</b>   | 50                 |

# GX Microscopes

product datasheet

## GXM Coverslips Specifications

### Coverslips



### Coverslips

| Type        | Dimensions mm | Thickness mm | Product Code | Qty per box |
|-------------|---------------|--------------|--------------|-------------|
| Square      | 18x18         | 0.13-0.17    | 18           | 100         |
| Square      | 20x20         | 0.13-0.17    | 20           | 100         |
| Square      | 22x22         | 0.13-0.17    | 22           | 100         |
| Square      | 24x24         | 0.13-0.17    | 24           | 100         |
| Rectangular | 24x32         | 0.13-0.17    | 24/32        | 100         |
| Rectangular | 24x40         | 0.13-0.17    | 24/40        | 100         |
| Rectangular | 24x50         | 0.13-0.17    | 24/50        | 100         |
| Rectangular | 24x60         | 0.13-0.17    | 24/60        | 100         |
| Circular    | 13            | 0.13-0.17    | 13C          | 100         |
| Circular    | 15            | 0.13-0.17    | 15C          | 100         |
| Circular    | 16            | 0.13-0.17    | 16C          | 100         |
| Circular    | 18            | 0.13-0.17    | 18C          | 100         |
| Circular    | 19            | 0.13-0.17    | 19C          | 100         |
| Circular    | 20            | 0.13-0.17    | 20C          | 100         |

Supplied to you by:

[www.gxmicroscopes.com](http://www.gxmicroscopes.com)