

GX Microscopes

product datasheet

V Series

Stands for Stereo
Microscopes




Specifications					
Model:	V1	V2	V3	V6	V7
Baseplate mm	254 x 254	Circular	G-Clamp	330 x 330	254 x 254
Main Pillar Diameter mm	38	38	38	32	38
Horizontal Pole Diameter mm	38	38	38	16	n/a
Total height mm	450	450	500	410	450
Horizontal Reach mm	276 - 573	276 - 573	130 - 390	370	240 - 846

The V series of stands for stereo microscopes provide a range of robust and convenient methods of mounting a wide range of stereo microscope heads.

The articulated V7 arm is particularly versatile with its long reach and 7 point articulation.

The stands are manufactured to a high quality and are designed for use in the most demanding environments.

Supplied to you by:



www.gxmicroscopes.com

Due to our policy of continuous development specifications may change without notice

June 04