

# GX Microscopes

## product datasheet

### Field Microscopes For Schools and Colleges



XSJ-20A



XTD-20A



XTD-40

The GX Microscopes range of college field microscopes are ideal for both field study and classroom work. There is a choice of models all at the same amazingly low price that will suit any budget.

The **XTD** models include adjustable binocular eyepieces and dual, angled objectives for stereoscopic 3D viewing. You can choose between the high magnification (x40) **XTD-40** model or the **XTD-20A** slightly lower magnification (x20) which benefits from a longer working distance (88mm).

For those students who prefer a simple monocular eyepiece the **XSJ-20A** is ideal, this model also benefits from an inclined eyepiece which can make prolonged observations more comfortable. The head can also be rotated to allow observation from the front or rear of the sample holder.

All of the microscopes are built to a robust, student-proof (as far as possible!) specification in metal, finished in stain and scratch resistant paint.

<b>Specifications:</b>			
<b>Model:</b>	<b>XSJ-20A</b>	<b>XTD-20A</b>	<b>XTD-40</b>
Total Magnification	20X	20X	40X
Objectives	2X	2X	4X
Eyepieces	WF10X		
Effective Working Distance	53mm	88mm	70mm
Head	Monocular Rotatable Inclined	Binocular adjustable	Binocular adjustable
Stand	Steel column pillar on a stable metal stand with slide holder clips and reversible background plate. Adjustable height.		

Supplied to you by:

[www.gxmicroscopes.com](http://www.gxmicroscopes.com)

Due to our policy of continuous development specifications may change without notice

June 04