

### UHM Universal Stands

UK manufactured long-arm stands for use with Meiji Stereo Microscopes, require F focussing ring/block and external illumination of specimens (see GX Microscope's illumination systems). Also suitable for use with other microscopes with the same microscope stand mountings.

**UHM** Universal stand with heavy base 8in x 10 in and SWIV 20mm swivel spigot mount to accept F focussing ring. Includes 40cm vertical pole with 40 cm horizontal arm.

**UHMFL** Universal stand with bolt-mounting table flange and SWIV 20mm swivel spigot mount to accept F focussing ring. Includes 40cm vertical pole with 40 cm horizontal arm.

**UHMSBU** Universal stand with 600mm deep x 400mm wide base-board and SWIV 20mm swivel spigot mount to accept F focussing ring. Includes 40cm vertical pole with 40 cm horizontal arm.

**UHMCB** Universal stand with heavy circular 9in base and counter balance arm



**UHM**